



# AVISCERA BIOSCIENCE

## Anti Rat Fetuin A Monoclonal IgG

### Product Information

Code	A00173-05-100
Name	Rat Fetuin A Mab
Clone No.	2B7E1B9D4
Lot No.	
Size	100 µg
Species	Rat
Host	Mouse
Immunogen	Rat Fetuin A rec.
Ab Type	IgG
Purification	Protein G
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Rat
Reconstitution	100 µl
Application	ELISA

AVISCERA BIOSCIENCE INC.  
2348 Walsh Ave. Suite C  
Santa Clara, CA 95051  
Tel: (408) 982 0300  
Fax: (408) 982 0301  
Email:  
Info@AvisceraBioscience.com  
www.AvisceraBioscience.com

### Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, E. coli-derived, recombinant Rat Fetuin A. That antibody was purified by Protein G affinity.

### Formulation

100 µg of Mouse anti Rat Fetuin A Monoclonal Antibody in 100 µl of PBS lyophilized form.

### Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (100µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted ( by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for 12 months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

### Specificity

This antibody has been selected for its ability to recognize recombinant Rat Fetuin A indirect ELISAs.

### Applications

**Indirect ELISA** - This antibody can be used at 1:8000 (0.125 µg/ml) to detect Rat Fetuin A on indirectly ELISA.

**ELISA Assay** - This antibody can be used as a capture antibody in a Rat Fetuin A sandwich immunoassay in combination with the Rat Fetuin A detection antibody (Code No.: A00173-02-100) and recombinant Rat Fetuin A (Code No.: 00173-02-100) as the standard. The suggested concentration range for this capture antibody is 2 µg/mL and should be titrated to determine the optimal concentration.

*Optimal dilutions should be determined by each laboratory for each application.*