



# AVISCERA BIOSCIENCE

## ANTI-HUMAN ANGIOPOIETIN LIKE 4 (ANGPTL4) IGG1

### Product Information

Code	A00309-01-50
Name	Human ANGPTL4 Mab
Clone No.	
Lot No.	
Size	50 µg
Species	Human
Host	Mice
Immunogen	ANGPTL4 (21- 250) (H), rec.
Ab Type	Mouse IgG1
Purification	Protein G
Formulation	PBS liquid
Carry	0.1% gelatin
Storage	-20 ° C
Specificity	Human
Application	ELISA

### Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with recombinant human ANGPTL (51-250). That antibody was purified by Protein G affinity.

### Formulation

50 µg of mouse anti human ANGPTL4 IgG1 in 250 µl of solution in phosphate-buffered saline (PBS) with 0.1% Gelatin .

### Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

### Specificity

This antibody has been selected for its ability to recognize human ANGPTL4 in western blot.

### Applications

**Western Blot** - This antibody can be used as antibody at 1-5 µg/ml with the appropriate detection antibody to detect human ANGPTL4 on transfected human tumor cells .

*Optimal dilutions should be determined by each laboratory for each application.*

AVISCERA BIOSCIENCE, INC  
2348 Walsh Ave, Suite C,  
Santa Clara, CA 95051  
USA  
Tel: (408) 982 0300  
Fax: (408) 982 0301  
info@Aviscerabioscience.com  
www.Aviscerabioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.