



AVISCERA BIOSCIENCE

Anti-Rat Adiponectin Antibody

Product Information

Code	A00010-06-100
Name	Rat Adiponectin Pab
Clone No.	Polyclonal
Lot No.	
Size	100 µg
Species	Rat
Host	Rabbit
Immunogen	Adiponectin, globular form (R), rec.
Ab Type	IgG
Purification	Protein A
Formulation	lyophilized Form without preservatives free
Carry	
Storage	-20 ° C
Specificity	Rat
Reconstitution	PBS, 100 µl
Application	IHC WB ELISA

AVISCERA BIOSCIENCE
2348 Walsh Ave., Suite c
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Info@Aviscerabioscience.com
www.AvisceraBioscience.com

Preparation

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant rat adiponectin, globular form. That was purified by Protein A affinity.

Formulation

100 µg of rabbit anti-rat Adiponectin antibody in PBS without preservatives was lyophilized.

Reconstitution

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize rat adiponectin in direct ELISAs and western blots as well as immunohistochemistry.

Applications

Direct ELISA - This antibody can be used at 1:5000 (0.20 µg/ml) with the appropriate secondary reagents to detect rat adiponectin, globular form.

Western blot - This antibody can be used at 1: 1000-500 [1-2 µg/mL] with the appropriate secondary reagents to detect adiponectin in rat adipose tissue at 30 KDa under reduce condition . The detection limit for rec. rat Adiponectin, globular from approximately at 10ng/lane under reducing conditions.

Immunohistochemistry-That Antibody can be used at 1: 500~1000 (1-2 µg/ml) with the appropriate secondary antibody to detect adiponectin in the paraffin embedded rat adipose tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.