



AVISCERA BIOSCIENCE

Rabbit Anti-Human Cerebellin-1 IgG

Product Information

Code	A00361-01-100
Name	CBLN-1 (H) Pab
Clone No.	N/A
Lot No.	
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	Cerebellin-1 (Human), mature form
Ab Type	IgG
Purification	Protein A
Formulation	Lyophilized Form without Preservatives
Carry	Free
Storage	-20 ° C
Specificity	H, R, M
Reconstitution	PBS, 100 µl
Application	IHC ELISA

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

Preparation

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant mature form of human Cerebellin-1. That IgG was purified by Protein A affinity.

Formulation

100µg of purified IgG in 100 µl of PBS without preservatives was lyophilized.

Reconstitution

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize mature form of human Cerebellin-1 on indirect ELISAs.

Applications

Indirect ELISA - This antibody can be used at 1: 6000 (0.166 µg/mL) with the appropriate secondary reagents to detect Cerebellin-1 (Human) on indirect ELISA.

Immunohistochemistry-That Antibody can be used at 2.5-5 µg /mL with the appropriate secondary antibody to detect Cerebellin-1 in human , mouse and rat brain tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.