



AVISCERA BIOSCIENCE

Rabbit Anti-Mouse CTRP12 Antibody

Product Information

Code	A00392-04-50B
Name	Anti- mouse CTRP12 Antibody
Clone No.	Polyclonal
Lot No.	
Size	50 µg
Species	Mouse
Host	Rabbit
Immunogen	Mouse CTRP12, GF rec.
Ab Type	IgG
Purification	Protein A
Formulation	Lyophilized form without preservatives free
Carry	
Storage	-20 ° C
Specificity	Mouse CTRP12
Reconstitution	50 µl
Application	IHC E

Preparation

This antibody was produced from a rabbit immunized with purified recombinant mouse CTRP12, Globular Form. That was purified by Protein A affinity.

Formulation

50 µg of rabbit anti-mouse CTRP12 IgG in 50 µl of TBS containing 0.2% BSA without preservatives in lyophilized form.

Reconstitution

Add 50 µl of PBS to the vial to prepare antibody stock solution at 50 µg/50 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Storage

This lyophilized antibody can be stored at 4 ° C for a few weeks. It can be stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Specificity

This antibody has been selected for its ability to recognize mouse CTRP12 on indirect ELISA and immunohistochemistry.

Applications

Indirect ELISA - This antibody can be used at 0.25 µg /mL with the appropriate secondary reagents to detect recombinant mouse CTRP12.

Immunohistochemistry-That Antibody can be used at 2-4 µg/ml with the appropriate secondary antibody to detect CTRP12 in mouse paraffin embedded visceral adipose tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Info@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.