



AVISCERA BIOSCIENCE

Anti Human CTRP5 (Complement C1q Tumor Necrosis Factor-Related Protein 5) Monoclonal Antibody, Biotinylated

Product Information

Code	A00594-07-50B
Name	Anti Human CTRP5 Monoclonal Antibody Biotinylated
Clone No.	7E11G5F9
Lot No.	
Size	50 µg
Species	Human
Host	Mouse
Immunogen	Human CTRP5 full length, rec.
Ab Type	Purified IgG
Purification	Protein G affinity
Formulation	lyophilized Form without preservatives
Carry	free
Storage	-20~ -70 ° C
Specificity	Human CTRP5
Reconstitution	PBS, 50 µl
Application	ELISA

ORDER INFORMATION

AVISCERA BIOSCIENCE, INC
2348 Walsh Ave., Suite C
Santa Clara, CA 95051

USA

Tel: (408) 982 0300

Info@Aviscerabioscience.com

www.Aviscerabioscience.com

PREPARATION

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified recombinant human CTRP5 full length. That antibody was purified Protein G affinity and conjugated with water soluble biotin.

FORMULATION

50 µg of Anti human CTRP5 Monoclonal antibody Biotinylated in 50 µl of PBS contains 0.01% BSA lyophilized form.

RECONSTITUTION

Add 50 µl of PBS to the vial to prepare antibody stock solution at 50 µg/50 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for up six months without detectable loss of activity.

STORAGE

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ~ -70 ° C for 12 months. **Avoid repeated freeze-thaw cycles.**

SPECIFICITY

This antibody has been selected for its ability to recognize human CTRP5 on indirect ELISA.

APPLICATIONS

Indirect ELISA - This antibody can be used at 1: 8000 (0.125 µg/ml) to detect recombinant Human CTRP5 on indirect ELISA.

Optimal dilutions should be determined by each laboratory for each application.

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.