



# AVISCERA BIOSCIENCE

## Mouse CTRP6 Globular Form Recombinant

### Product Information

Code	00089-04-100
Name	Mouse CTRP6 Globular Form, rec.
Lot No.	
Size	100 µg
Species	Mouse
Sequence	Globular form
Protein ID	Q6IR41
Gene ID	NM_028331
MW	32 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel PBS lyophilized
Formulation	Form without preservatives
Carry	free
Storage	-20 ° C ~ -70 ° C
Reconstitution	500 µl
Application	ELISA Cell biology

### Description

A DNA sequence encoding the mouse CTRP6 , globular form (Cys<sup>125</sup>-Asn<sup>264</sup>) with 6 His tag on the N-Terminus was expressed in *E. Coli* . This protein was purified by Ni-NTA column.

### Formulation

Lyophilized 100 µg mouse CTRP6 globular form in 100 µl of PBS, Carry free.

### Reconstitution & Storage

Add 500 µl deionized water to the vial to prepare a working stock solution at 200 µg/mL . Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at - 20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

**Sequence:** mouse CTRP6 , globular form (Cys<sup>125</sup>-Asn<sup>264</sup>)

```
25 CQTHYS AFSVGRKTGL HSENFLSLL FDRVFNVDG
HFDMATGSFV APLRGLYFFS LNVHSWNYKE TYVHIVHNEQ
AVVILYAQPS ERSIMQSQSV MLPLVPGDRV WVRLFKRERE
NGIYSDDVDVT YITFSGHLLIK AEDN 264
```

AVISCERA BIOSCIENCE

2348 Walsh Ave, Suite C  
Santa Clara, CA 95051  
USA

Tel: (408) 982 0300

Fax: (408) 982 0301

Info@Aviscerabioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.