



# AVISCERA BIOSCIENCE

## Rabbit Anti- Human CTRP6 IgG

### Product Information

Code	A00089-01-100
Name	Human CTRP6 Pab
Clone No.	N/A
Lot No.	
Size	100 µl
Species	Human
Host	Rabbit
Immunogen	CTRP6 (Human) mature form rec.
Ab Type	IgG
Purification	Protein A Lyophilized
Formulation	Form without Preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	PBS, 100 µl
Application	IHC ELISA

### Preparation

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant mature form of human CTRP6. That IgG was purified by Protein A affinity

### Formulation

100µg of purified IgG in 100 µl of PBS without preservatives was lyophilized.

### Reconstitution

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted ( by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

### Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

### Specificity

This antibody has been selected for its ability to recognize mature form human REG-1α in direct ELISAs as well as immunohistochemistry.

### Applications

**Indirect ELISA** - This antibody can be used at 1: 4000 (0.25 µg /mL) with the appropriate secondary reagents to detect human CTRP6 mature form on indirect ELISA.

**Immunohistochemistry**-That Antibody can be used at 1: 200 (5 µg /mL) with the appropriate secondary antibody to detect CTRP6 in human adipose and heart tissues (ABC).

*Optimal dilutions should be determined by each laboratory for each application.*

AVISCERA BIOSCIENCE  
2348 Walsh Ave., Suite C  
Santa Clara, CA 95051  
USA  
Tel: (408) 982 0300  
Info@AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.