



# AVISCERA BIOSCIENCE

## Anti-Rat Fetuin A IgG Biotinylated

### Product Information

Code	A00173-02-50B
Name	Anti Rat Fetuin A IgG Biotinylated
Clone No.	N/A
Lot No.	
Size	50 µg
Species	Rat
Host	Rabbit
Immunogen	Rat Fetuin A rec.
Ab Type	IgG
Purification	Protein A affinity
Formulation	lyophilized Form without preservatives free
Carry	
Storage	-20~ -70 ° C
Specificity	Rat
Reconstitution	PBS, 50 µl
Application	ELISA

### PREPARATION

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant rat Fetuin A. That IgG was purified by Protein A affinity and conjugated with water soluble biotin.

### FORMULATION

50 µg of Rabbit Anti Rat Fetuin A purified IgG Biotinylated in 50 µl of PBS contains 0.1% BSA without preservatives was lyophilized.

### RECONSTITUTION

Add 50 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for up six months without detectable loss of activity.

### STORAGE

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ~ -70 ° C for 12 months. **Avoid repeated freeze-thaw cycles.**

### SPECIFICITY

This antibody has been selected for its ability to recognize rat Fetuin A in indirect ELISA.

### APPLICATIONS

**Indirect ELISA** - This antibody can be used at 1:10000 (0.1µg/ml) to detect recombinant rat Fetuin A on indirect ELISA.

**ELISA**-This antibody can be used as detection antibody at 0.25-0.5 µg/ml combines with a capture antibody A00173-02-100 to detect rat fetuin A recombinant 00173-02-100 on ELISA.

### ORDER INFORMATION

AVISCERA BIOSCIENCE, INC  
2348 Walsh Ave., Suite C  
Santa Clara, CA 95051  
USA  
Tel: (408) 982 0300  
Info@Aviscerabioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.