



AVISCERA BIOSCIENCE

Anti Rat Fetuin A Monoclonal IgG

Product Information

Code	A00173-08-100
Name	Rat Fetuin A Mab
Clone No.	2F5H3G3C9
Lot No.	
Size	100 µg
Species	Rat
Host	Mouse
Immunogen	Rat Fetuin A rec.
Ab Type	IgG
Purification	Protein G
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Rat
Reconstitution	100 µl
Application	ELISA

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, E. coli-derived, recombinant Rat Fetuin A. That antibody was purified by Protein G affinity.

Formulation

100 µg of Mouse anti Rat Fetuin A Monoclonal Antibody in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (100µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for 12 months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize recombinant Rat Fetuin A indirect ELISAs.

Applications

Indirect ELISA - This antibody can be used at 1:1000 (1 µg/ml) to detect Rat Fetuin A on indirectly ELISA.

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Info@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.