



AVISCERA BIOSCIENCE

Human Hepatoma-Derived Growth Factor, Related Protein 3 (HDGFRP3), Recombinant

Product Information

Code	00305-01-100
Name	HDGFRP3, rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	Met1-Thr203
Protein ID	Q9Y3E1
Gene ID	51810
MW	22.6 KDa
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel PBS lyophilized
Formulation	Form without preservatives
Carry	Free
Storage	-20 ° C ~ -70° C
Reconstitution	500 µl
Application	ELISA

Description

A DNA sequence encoding the human HDGFRP3 (Met¹-Thr²⁰³) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 100 µg human HDGFRP3 (Met¹-Thr²⁰³) in 100 µl of PBS. Carry free.

Reconstitution & Storage

Add 500 µl deionized water to the vial to prepare a working stock solution at 200 µg/mL. Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: Human HDGFRP3 (Met¹-Thr²⁰³)

```
1
MARPRPREYK AGDLVFAKMK GYPHWPARID ELPEGAVKPP
ANKYPIFFFG THETAFLGPK DLFPYKEYKD KFGKSNKRKG
FNEGLWEIEN NPGVKFTGYQ AIQQQSSSET EGEGGNTADA
SSEEEGDRVE EDGKGRKNE KAGSKRKKSY TSKKSSKQSR
KSPGDEDDKD CKEEENKSSS EGGDAGNDTR NTTSDLQKTS
EGT 203
```

AVISCERA BIOSCIENCE
2348 Walsh Ave., Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Info@AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.