



AVISCERA BIOSCIENCE

Mouse Meteorin Like Protein (METRNL) Recombinant

Product Information

Code	00478-03-100
Name	Meteorin Like Protein (Mouse), rec.
Lot No.	
Size	100 µg
Species	Mouse
Sequence	Mature form
Protein ID	NP_659046
Gene ID	NM_144797
MW	32 KDa
Source	E. Coli
Purity	>95% in SDS- PAGE gel PBS lyophilized
Formulation	Form without preservatives
Carry	free
Storage	-20 °C ~ -70° C
Reconstitution	500 µl
Application	Cell Biology

Alternative name(s):

Subfatin

Description

A DNA sequence encoding the mature mouse Meteorin Like Protein (Gln⁴⁶-Glu³¹¹) with 6 x His Tag on the N-Terminus was expressed in *E. Coli*.

Formulation

Lyophilized 100 µg mouse Meteorin Like Protein in 100 µl PBS, pH7.4, carry free.

Reconstitution & Storage

Reconstitute at 500 µl (200 µg/mL) in sterile PBS containing at least 0.1% BSA. Allow to set at least 30 minutes at 4 °C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 12 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: mature form mouse Meteorin Like Protein (Q46-E311)

ORDER INFORMATION
AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Info@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.