



AVISCERA BIOSCIENCE

ANTI-HUMAN EG-VEGF/PROKINETICIN 1 IGG BIOTINYLATED

Product Information

Code	A00227-01-50B
Name	Human EG-VEGF IgG Biotinylated
Clone No.	Polyclonal
Lot No.	
Size	50 µg
Species	Human
Host	Rabbit
Immunogen	Human EG-VEGF , rec.
Ab Type	IgG
Purification	Protein A
Formulation	lyophilized Form without preservatives free
Carry	
Storage	-20 ° C
Specificity	Human
Reconstitution	50 µl
Application	IHC ELISA

Preparation

This antibody was produced from a rabbit immunized with purified recombinant human EG-VEGF. That was purified by Protein A affinity and conjugated with water soluble biotin.

Formulation

50 µg of rabbit anti-human EG-VEGF IgG in 50 µl of TBS containing 0.2% BSA without preservatives was lyophilized.

Reconstitution

Add 50 µl of TBS to the vial to prepare antibody stock solution at 50 µg/50 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize human EG-VEGF in direct ELISAs and immunohistochemistry.

Applications

Indirect ELISA - This antibody can be used at 0.25 µg/ml with the appropriate secondary reagents to detect human EG-VEGF.

Immunohistochemistry-That Antibody can be used at 2-4 µg/ml with the appropriate secondary antibody to detect EG-VEGF in the paraffin embedded human colon cancer tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE, INC
2348 WALSH AVE., SUITE C
SANTA CLARA, CA 95051
USA
TEL: (408) 982 0300
Info@Aviscerabioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

AVISCERA BIOSCIENCE