

HUMAN VISCERAL ADIPOSE TISSUE-DERIVED SERINE PROTEASE INHIBITOR (VASPIN) ELISA KIT

FOR THE QUANTITATIVE DETERMINATION OF
HUMAN VASPIN CONCENTRATIONS IN SERUM
AND PLASMA



FOR RESEARCH USE ONLY. NOT FOR USE IN
DIAGNOSTIC PROCEDURES.

PURCHASE INFORMATION:

ELISA NAME	HUMAN VASPIN ELISA
Catalog No.	SK00560-01
Lot No.	
Formulation	96 T
Standard range	0.312-20 ng/ml
Sensitivity	0.156 ng/ml
Sample Volume	100 μ l
Sample Type	Serum, EDTA Plasma
Dilution Factor	Optimal dilutions should be determined by each laboratory for each application
Specificity	Human Vaspin
Intra-assay Precision	4-6%
Inter-assay Precision	8-10%
Storage	4 °C

ORDER CONTACT:
AVISCIERA BIOSCIENCE INC.
 2348 Walsh Ave., Suite C
 Santa Clara, CA 95051
 Tel: (408) 982 0300
 Fax: (408) 982 0301
 Email: Sales@AvisceraBioscience.com
Info@AvisceraBioscience.com
www.AvisceraBioscience.com

INTRODUCTION

Human VASPIN immunoassay is a 3.5 - 4.5 hour solid phase ELISA designed to measure human VASPIN in serum and plasma. It contains recombinant human VASPIN and antibodies raised against this protein. It has been shown to accurately quantify recombinant human VASPIN. Results obtained with naturally occurring VASPIN samples showed linear curves that were parallel to the standard curves obtained using the kit standards. These results indicate that the immunoassay kit can be used to determine relative mass values for natural human VASPIN.

PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for VASPIN has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any VASPIN present is bound by the immobilized antibody. After washing away any unbound substances, an antibody specific for VASPIN is added to the wells. Following a wash to remove any unbound detection antibody reagent, an anti-rabbit IgG HRP conjugate is added to the wells. After washing away any unbound enzyme, a substrate solution is added to the wells and color develops in proportion to the amount of VASPIN bound in the initial step. The color development is stopped and the intensity of the color is measured.

LIMITATIONS OF THE PROCEDURE

_ FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

_ The kit should not be used beyond the expiration date on the kit label.

_ Do not mix or substitute reagents with those from other lots or sources.

_ It is important that the Dilution Buffer selected for the standard curve be consistent with the samples being assayed.

_ If samples generate values higher than the highest standard, dilute the samples with Dilution Buffer and repeat the assay.

_ Any variation in standard diluent, operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

_ This assay is designed to eliminate interference by soluble receptors, binding proteins, and other factors present in biological samples. Until all factors have been tested in the immunoassay, the possibility of interference cannot be excluded.

MATERIALS PROVIDED

DESCRIPTION	CODE	QUANTITY
VASPIN Microplate - 96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against VASPIN.	560-01-01	1 plate
VASPIN Standard – 10 ng/vial of recombinant human VASPIN in a buffered protein base with preservatives; lyophilized.	560-01-02	2 vials
Detection Antibody Concentrate – 105 µL/vial, 100-fold concentrated of an antibody against VASPIN with preservatives; lyophilized.	560-01-03	1 vial
Positive Control - one vial of recombinant human VASPIN, lyophilized	560-01-04	1 vial
Anti-Rabbit IgG HRP Conjugate - 120 µl/vial, 100-fold concentrated solution of Anti Rabbit IgG HRP conjugate to HRP	ARIGHRP	1 vial
Dilution Buffer - 60mL of buffered protein based solution with preservatives	DB08	1 bottle
Wash Buffer -50 mL of 10-fold concentrated buffered surfactant, with preservative.	WB01	1 bottle
TMB Substrate Solution - 11 mL of TMB substrate solution	TMB01	1 bottle
Stop Solution – 11 mL of 0.5M HCl solution	S-STOP	1 bottle
Plate Sealer	EAPS	1 piece
Plastic Pouch	P01	1 piece

STORAGE

Unopened Kit: Store at 2 - 8° C for up to 6 months. For longer storage, unopened Standard, Positive Control and Detection Antibody Concentrate should be stored at -20°C or -70°C. Do not use kit past expiration date.

Opened / Reconstituted Reagents: Reconstituted Standard and Antibody Solution SHOULD BE STORED at -20°C or -70°C for up to one month. ARIG-HRP Conjugate 100-fold concentrated and other components may be stored at 2 - 8°C for up to 6 months.

Microplate Wells: Return unused wells to the plastic pouch with the desiccant pack and seal along the entire edge of zip-seal. Microplate may be stored for up to 6 months at 2 - 8°C.

OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Microplate shaker (250-300rpm).
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 100 mL and 500 mL graduated cylinders.

SAMPLE COLLECTION AND STORAGE

Plasma - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤-20°C. Avoid repeated freeze-thaw cycles.

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at ≤ -20°C. Avoid repeated freeze-thaw cycles.

Note: Use Aprotinin (enzyme inhibitor) (Code No.00700-01-25) for ALL sample collection to prevent sample degradation. 0.5 TIU per ml of sample solution.

REAGENT PREPARATION


Bring all reagents to room temperature before use.

Wash Buffer - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 50 mL of Wash Buffer Concentrate into deionized or distilled water (450 mL) to prepare 500 mL of Wash Buffer.

VASPIN Standard - Refer to vial label for reconstitution volume. Reconstitute the VASPIN standard with 0.5 mL of Dilution Buffer. This reconstitution produces a stock solution of 20 ng/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Pipette 250 µL of Dilution Buffer into tubes #1 to #6. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 20 ng/mL standard serves as the

high standard. The Dilution Buffer serves as the zero standard (0 ng/mL).

TUBE	STANDARD	DILUTION BUFFER	CONCENTRATION
stock	Powder	0.5 ml	20 ng/ml
# 1	250µl of stock	250µl	10 ng/ml
# 2	250µl of 1	250µl	5 ng/ml
# 3	250µl of 2	250µl	2.5 ng/ml
# 4	250µl of 3	250µl	1.25 ng/ml
# 5	250µl of 4	250µl	0.625 ng/ml
# 6	250µl of 5	250µl	0.3125 ng/ml



Detection Antibody Concentrate - Reconstitute the Detection Antibody Concentrate with 105 µL of Dilution Buffer to produce a 100-fold concentrated stock solution. Pipette 10.395 mL of Dilution Buffer into a 15 mL centrifuge tube and transfer 105 µL of 100-fold concentrated stock solution to prepare working solution.

Anti-Rabbit IgG HRP Conjugate - Pipette 11.88 mL of Dilution Buffer into a 15 mL centrifuge tube and transfer 120 µL of 100-fold concentrated stock solution to prepare working solution. **Note:** 1x working solution of Anti-Rabbit IgG HRP Conjugate should be used within a few days.

Positive Control - Reconstitute the Positive Control with 1.0 mL of Dilution Buffer. **Note:** Positive Control should be prepared and used within a few days.

ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that

blank, standards, positive control and samples be assayed in duplicates.


1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess micro-plate strips from the plate frame, return them to the plastic pouch with the desiccant pack and seal.
3. Add 100 µL of Dilution Buffer to Blank wells (B4, B5).
4. Add 100 µL of Standard (from B2, B3 to H2, H3), sample, or positive control (D4, D5) per well. Cover with plate sealer. Incubate for 2 hours on micro-plate shaker at room temperature. A plate layout is provided to record standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with 1x Wash Buffer (300 µL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100 µL of Detection Antibody working solution to each well. Cover with plate sealer. Incubate for 2 hours on micro-plate shaker at room temperature.
7. Repeat the aspiration/wash as in step 5.
8. Add 100 µL of Anti-Rabbit IgG HRP Conjugate working solution to each well. Incubate for 1 hour on micro-plate shaker at room temperature.
Protect from light.
9. Repeat the aspiration/wash as in step 5.
10. Add 100 µL of Substrate Solution to each well. Incubate for 8-12 minutes at room temperature.
Protect from light.
11. Add 100 µL of Stop Solution to each well. The color in the wells should change from blue to yellow. If the color in the wells is green, or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
12. Determine the optical density of each well within 15 minutes, using a micro-plate reader set to 450 nm.

CALCULATION OF RESULTS

Average the duplicate readings for each standard, positive control, and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a log-log curve fit. As

an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the VASPIN concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.



TYPICAL DATA

This standard curve data is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

STANDARD (NG/ML)	AVERAGE OD450 (CORRECTED)
Blank	0 (0.064)
0.156 (optional)	0.004
0.3125	0.014
0.625	0.042
1.25	0.078
2.5	0.157
5	0.638
10	1.172
20	1.707

- Lot No.:
- Positive Control: 2.3-4.3 ng/mL

CALIBRATION

This immunoassay is calibrated against a highly purified recombinant human VASPIN.

SENSITIVITY

Twenty-five assays were evaluated and the minimum detectable dose (MDD) of VASPIN was 0.156 ng/mL.

SPECIFICITY

This assay recognizes both natural and recombinant human Vaspin. The factors listed below were prepared at 200 ng/mL in Dilution Buffer, and assayed for cross reactivity. No significant cross-reactivity or interference was observed.

PROTEINS	CROSS-REACTIVITY
Human VASPIN	100%
Mouse VASPIN	0
Human Omentin 1	0
Human Visfatin	0
Human Leptin	0
Human CTRP1	0
Human CTRP9	0
Human SPARC	0

LINEARITY

To assess the linearity of the assay , pooled research serum samples containing human Vaspin was diluted in Dilution Buffer and then assayed.

DILUTION FACOR	ASSAYED (NG/ML)	FINAL (NG/ML)	RECOVERY (%)
1 x	2.072	2.072	100
2 x	1.102	2.204	106.4
4 x	0.506	2.024	97.7

SUMMARY OF ASSAY PROCEDURE

