



AVISCERA BIOSCIENCE

Mouse Betatrophin /Angiopoietin-like Protein 8 /Lipasin Recombinant

Product Information

Code	00528-03-100
Name	Mouse Betatrophin
Lot No.	
Size	100 µg
Species	Mouse
Sequence	Mature form
Protein ID	NP_001074409
Gene ID	NM_001080940
MW	25 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel
Formulation	PBS lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C ~ -70° C
Reconstitution	100 µl
Application	ELISA Cell biology

Description

A DNA sequence encoding the mature form of mouse Betatrophin (Val¹⁶-Ala¹⁹⁸) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 100 µg mouse Betatrophin in 100 µl of PBS, Carry free.

Reconstitution & Storage

Add 100 µl PBS to the vial to prepare a working stock solution at 1000 µg/mL . Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at - 20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: mature form of mouse Betatrophin (Val¹⁶-Ala¹⁹⁸)

AVISCERA BIOSCIENCE

2348 Walsh Ave, Suite C

Santa Clara, CA 95051

USA

Tel: (408) 982 0300

Fax: (408) 982 0301

Info@AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.