



AVISCERA BIOSCIENCE

Human Betatrophin/Hepatocellular Carcinoma-associated Gene TD26 Isoform CRA_a Recombinant

Product Information

Code	00528-07-100
Name	Human Betatrophin Isoform CRA_a
Lot No.	
Size	100 µg
Species	Human
Sequence	M1-A99
Protein ID	EAW84185
Gene ID	55908
MW	14.43 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel
Formulation	PBS lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C ~ -70° C
Reconstitution	100 µl
Application	ELISA Cell biology

Description

A DNA sequence encoding the mature form of human Betatrophin /HCAG TD26 Isoform CRA_a (Met¹-Ala⁹⁹) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 100 µg human Betatrophin Isoform CRA_a recombinant in 100 µl of PBS, Carry free.

Reconstitution & Storage

Add 100 µl TBS, PH8.0 to the vial to prepare a working stock solution at 1000 µg/mL . Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at - 20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: human Betatrophin Isoform CRA_a (Met¹-Ala⁹⁹)

AVISCERA BIOSCIENCE

2348 Walsh Ave, Suite C

Santa Clara, CA 95051

USA

Tel: (408) 982 0300

Fax: (408) 982 0301

Info@AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.