



AVISCERA BIOSCIENCE

Human Angiopoietin Like 3 (ANGPTL3) Recombinant

Product Information

Code	00311-01-100
Name	Human ANGPTL3 Rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	Ser17-Glu460
Protein ID	NP_055310
Gene ID	NM_014495
MW	55.6 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS- PAGE gel PBS lyophilized form without preservatives
Formulation	
Carry	Free
Storage	-20 ° C ~ -70° C
Reconstitution	100 µl
Application	Cell biology

Description

A DNA sequence encoding the mature form of human ANGPTL3 (Ser¹⁷-Glu⁴⁶⁰) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 100 µg of human Angiopoietin Like 3 in 100 µl of PBS Carry free.

Reconstitution & Storage

Add 100 µl deionized water to the vial to prepare a working stock solution at 1000 µg/mL. Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: Human Angiopoietin Like 3 (Ser¹⁷-Glu⁴⁶⁰)

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.