



# AVISCERA BIOSCIENCE

## Anti Albumin (Human) Monoclonal Antibody (Capture)

### Product Information

Code	A00383-07-100
Name	Anti calreticulin Mab
Clone No.	2B9
Lot No.	
Size	100 µg
Species	Human
Host	Mouse
Immunogen	Human albumin.
Ab Type	IgG
Purification	Protein G
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	100 µl
Application	ELISA

### Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, human albumin. That antibody was purified by Protein G affinity.

### Formulation

100 µg of Anti albumin (Human) Monoclonal Antibody in 100 µl of PBS lyophilized form.

### Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (100µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted ( by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for 12 months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

### Specificity

This antibody has been selected for its ability to recognize albumin (Human) on ELISA.

### Applications

**ELISA**-This antibody can be used as capture antibody at 1-2 µg/ml combine the detection antibody A00383-01-100B to detect human albumin on ELISA.

*Optimal dilutions should be determined by each laboratory for each application.*

AVISCERA BIOSCIENCE INC.  
2348 Walsh Ave. Suite C  
Santa Clara, CA 95051  
USA  
Tel: (408) 982 0300  
Fax: (408) 982 0301  
Email:  
Sales@AvisceraBioscience.com  
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.