



AVISCERA BIOSCIENCE

Rabbit Anti-Human BMP-5 IgG

Product Information

Code	A00013-01-100
Name	Human BMP5 Pab
Clone No.	N/A
Lot No.	HBMP5001
Size	100 µl
Species	Human
Host	Rabbit
Immunogen	BMP-5(H) rec.
Ab Type	IgG
Purification	Protein A
Formulation	lyophilized Form without preservatives free
Carry	
Storage	-20 ° C
Specificity	Human
Reconstitution	PBS, 100 µl
Application	ELISA

Preparation

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant Human BMP-5, mature form. That IgG was purified by Protein A affinity

Formulation

100 µg of purified IgG in PBS without preservatives was lyophilized.

Reconstitution

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 uL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize recombinant human BMP-5 in direct ELISA

Applications

Direct ELISA - This antibody can be used at 1: 5000 (0.2 µg /mL) with the appropriate secondary reagents to detect human BMP-5.

Optimal dilutions should be determined by each laboratory for each application.

ORDER INFORMATION

AVISCERA BIOSCIENCE, INC
2348 Walsh Ave., Suite C
Santa Clara, CA 95051
USA
Info@Aviscerabioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.