



Human Cerebellin-1 (CBLN1), Recombinant

Product Information

Code	00361-01-100
Name	Cerebellin-1 (H), rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	mature
Protein ID	P23435
Gene ID	869
MW	33 KD
Tag	His tag on N terminal
Source	<i>E. Coli</i>
Purity	>95% in SDS-PAGE gel
Carry	free
Storage	-20 ° C ~ -70° C
Application	Cell Biology

Description

A DNA sequence encoding the human Cerebellin-1 (Gln²²-Leu¹⁹³) with 6 His tag on the N-Terminus was expressed in *E. Coli*. and purified by Ni-NTA column

Formulation

Lyophilized 100 µg human Cerebellin-1 from a 100 µl of PBS. Carry free.

Reconstitution & Storage

Add 500 µl sterile deionized water to the vial to prepare a working stock solution at 200 µg/mL . Allow to set at least 30 minutes at 4 ° C, mix well.

Lyophilized protein is stable for up to 8 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: Extracellular domain of human Cerebellin-1 (Gln²²-Leu¹⁹³)

```
22
  QNETEPIVL  EGKCLVVCDS  NPTSDPTGTA  LGISVRSRSA
KVAFSAIRST  NHEPSEMSNR  TMIIYFDQVL  VNIGNNFDSE
RSTFIAPRKG  IYSFNFHVVK  VYNRQTIQVS  LMLNGWPVIS
AFAGDQDVTR  EAASNGVLIQ  MEKGDRAVYK  LERGNLMGGW
KYSTFSGFLV  FPL 193
```

AVISCERA BIOSCIENCE, INC
2348 WALSH AVE., SUITE C
SANTA CLARA, CA 95051
USA

Tel: (408) 982 0300

Fax: (408) 982 0301

Info@AvisceraBioscience.com

www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.