



Anti CTRP13 (Human, Mouse) Monoclonal Antibody

Product Information

Code	A00333-04-100
Name	CTRP13 (Human, Mouse) Mab
Clone No.	1D8H6B4
Lot No.	
Size	100 µg
Species	Human
Host	Mouse
Immunogen	Human ctrp13 rec.
Ab Type	IgG
Purification	Protein G
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human, Mouse
Reconstitution	100 µl
Application	ELISA

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant ctrp13 (Human, Mouse) mature form. That antibody was purified by Protein G affinity.

Formulation

100 µg of Anti CTRP13 (Human) Monoclonal Antibody in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (100µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for 12 months without detectable loss of activity.

Avoid repeated freeze-thaw cycles.

Specificity

This antibody has been selected for its ability to recognize recombinant CTRP13 (Human, Mouse) in indirect ELISAs.

Applications

Indirect ELISA - This antibody can be used at 1:25600 (39 ng/ml) to detect CTRP13 (Human, Mouse) on indirectly ELISA.

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com