



AVISCERA BIOSCIENCE

Human CTRP7/C1QTNF7 Globular Form Recombinant

Product Information

Code	00396-02-100
Name	CTR7 Globular Form (Human), rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	S143-L289
Protein ID	Q9BXJ2
Gene ID	114905
MW	16 KD
Tag	His tag on N terminus
Source	E. Coli
Purity	>90% in SDS- PAGE gel PBS
Formulation	lyophilized form without preservatives
Carry	free
Storage	-20 °C ~ -70 °C
Reconstitution	100 µl
Application	ELISA Cell biology

C1QTNF7/Complement C1q tumor necrosis factor-related protein 7

Description

A DNA sequence encoding the human CTRP7, Globular Form (Ser¹⁴³-Leu²⁸⁹) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 100 µg human CTRP7 Globular Form in PBS. Carry free.

Reconstitution & Storage

Add 100 µl PBS or deionized water to the vial to prepare a working stock solution at 100 µg/100 µl. Allow to set at least 30 minutes at 4 °C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: human CTRP7, Globular Form (Ser¹⁴³-Leu²⁸⁹)

143

```
SIVLKSAF SVGITTSYPE ERLPIIFNKV LFNEGEHYNP  
ATGKFICAFP GIYYFSYDIT LANKHLAIGL VHNGQYRIKT  
FDANTGNHDV ASGSTVIYLLQ PEDEVWLEIF FTDQNGLFSD  
PGWADSLFSG FLLYVDTDYL DSISEDDEL 289
```

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

ORDER INFORMATION

Aviscera Bioscience Inc.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Info@AvisceraBioscience.com