



AVISCERA BIOSCIENCE

Rabbit Anti-Human CTRP7 IgG Biotinylated

Product Information

Code	A00396-02-50B
Name	Anti- Human CTRP7 IgG Biotinylated
Clone No.	Polyclonal
Lot No.	
Size	50 µg
Species	Human
Host	Rabbit
Immunogen	H CTRP7, GF, rec.
Ab Type	IgG
Purification	Protein A
Formulation	Lyophilized form without preservatives free
Carry	
Storage	-20 ° C
Specificity	Human CTRP7
Reconstitution	50 µl
Application	IHC E

Preparation

This antibody was produced from a rabbit immunized with purified recombinant human CTRP7, Globular form. That was purified by Protein A affinity and conjugated with water soluble biotin.

Formulation

50 µg of rabbit anti-Human CTRP7 IgG Biotinylated in 50 µl of TBS contains 0.2% BSA without preservatives in lyophilized form.

Reconstitution

Add 50 µl of TBS to the vial to prepare antibody stock solution at 50 µg/50 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Storage

This lyophilized antibody can be stored at 4 ° C for a few weeks. It can be stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Specificity

This antibody has been selected for its ability to recognize human CTRP7 on indirect ELISA and immunohistochemistry.

Applications

Indirect ELISA - This antibody can be used at 0.125 µg /mL with the appropriate secondary reagents to detect recombinant human CTRP7.

Immunohistochemistry-That Antibody can be used at 2-4 µg/ml with the appropriate secondary antibody to detect CTRP7 in human paraffin embedded visceral adipose tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Info@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.