



## Anti Human CTRP7 (Complement C1q Tumor necrosis factor-Related Protein 7) Monoclonal Antibody Biotinylated

### Product Information

Code	A00396-06-50B
Name	Human CTRP7 Mab Biotinylated
Clone No.	2C7F6F4
Lot No.	
Size	50 µg
Species	Human
Host	Mouse
Immunogen	Human gCTRP7, rec.
Ab Type	Purified IgG
Purification	Protein G affinity
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	100 µl
Application	ELISA

AVISCERA BIOSCIENCE INC.  
2348 Walsh Ave. Suite C  
Santa Clara, CA 95051  
Tel: (408) 982 0300  
Fax: (408) 982 0301  
Email:  
Sales@AvisceraBioscience.com  
www.AvisceraBioscience.com

### Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified recombinant human CTRP7, globular form, His Tag on N-Terminal. That antibody was purified Protein G affinity and conjugated with water soluble biotin.

### Formulation

100 µg of Anti human CTRP7 monoclonal antibody biotinylated in 100 µl of PBS lyophilized form.

### Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (500ug/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted ( by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for six months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

### Specificity

This antibody has been selected for its ability to recognize human CTRP7 globular form in indirect ELISAs. But it does not show any crossreactivity with human gCTRP3, gCTRP9, gCTRP4, CTRP2, CTRP10 and gCTRP6.

### Applications

**Indirect ELISA** - This antibody can be used at 0.25 µg/ml to detect human gCTRP1 on indirectly ELISA.

*Optimal dilutions should be determined by each laboratory for each application.*