



AVISCERA BIOSCIENCE

Human Cyclophilin A (CYPA) Recombinant

Product Information

Code	00558-01-100
Name	Human CYPA,rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	V2-E165
Protein ID	NP_066953.1
Gene ID	NM_021130.4
MW	21 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS- PAGE gel PBS lyophilized
Formulation	Form without preservatives
Carry	free
Storage	-20 ° C ~ -70° C
Reconstituti on	500 µl

Description

A DNA sequence encoding the mature form of human Cyclophilin A (V²-E¹⁶⁵) with 6 His tag on the N-Terminus was expressed in *E. Coli* . This protein was purified by Ni-NTA column.

Formulation

Lyophilized 100 µg human Cyclophilin A in 100 µl of PBS, Carry free.

Reconstitution & Storage

Add 500 µl deionized water to the vial to prepare a working stock solution at 200 µg/mL . Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at - 20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: Human Cyclophilin A (V2-E165)

AVISCERA BIOSCIENCE

2348 Walsh Ave, Suite C
Santa Clara, CA 95051
USA

Tel: (408) 982 0300

Fax: (408) 982 0301

Info@Aviscerabioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.