



Human Fatty Acid Binding Protein-7 (FABP-7/BFABP) Monoclonal Antibody

Product Information

Code	A00191-03-100
Name	Human FABP-7 MAb
Clone No.	1G1B5
Lot No.	
Size	100 µg
Species	Human
Host	Mouse
Immunogen	Human FABP-7, rec.
Ab Type	IgG
Purification	Protein G
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 °C
Specificity	Human
Reconstitution	N/A
Application	ELISA IHC

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, E. coli-derived, recombinant human FABP-7. That antibody was purified by Protein G affinity.

Formulation

100 µg of Mouse IgG in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (100µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 uL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for six months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize human FABP-7 on indirect ELISA.

Applications

Indirect ELISA - This antibody can be used at 0.1 µg/ml to detect recombinant human FABP-7 on indirectly ELISA.

Immunohistochemistry- This antibody can be used at 1:50 (20 µg/ml) with the appropriate secondary antibody to detect FABP-7 in human brain tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Info@AvisceraBioscience.com
www.AvisceraBioscience.com