



AVISCERA BIOSCIENCE

Human Pancreatic-derived Factor (FAM3B), Recombinant

Product Information

Code	00344-01-100
Name	FAM3B (H), rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	mature
Protein ID	P58499
Gene ID	54097
MW	40 KD
Tag	His tag on N terminal
Source	<i>E. Coli</i>
Purity	>95% in SDS-PAGE gel
Carry	free
Storage	-20 ° C ~ -70° C
Application	Cell Biology

Description

A DNA sequence encoding the human Pancreatic-derived Factor (Glu³⁰-Ser²³⁵) with 6 His tag on the N-Terminus was expressed in *E. Coli*. and purified by Ni-NTA column.

Formulation

Lyophilized 100 µg human Pancreatic-derived factor (FAM3B) from a 50 µl Tris Buffer solution. Carry free.

Reconstitution & Storage

Add 500 µl TBS to the vial to prepare a working stock solution at 200 µg/mL . Allow to set at least 30 minutes at 4 ° C, mix well. Lyophilized protein is stable for up to 8 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: human Pancreatic-derived Factor (Glu³⁰-Ser²³⁵)

```

                                                                 30E
LIPDAPLSSA  AYSIRSIGER  PVLKAPVPKR  QKCDHWTPCP
SDTYAYRLLS  GGGRSKYAKI  CFEDNLLMGE  QLGNVARGIN
IAIVNYVTGN  VTATRCFDY  EGDNSGPMTK  FIQSAAPKSL
LFMVTYDDGS  TRLNNDAKNA  IEALGSKEIR  NMKFRSSWVF
IAAKGLELPS  EIQRKINHS  DAKNNRYSGW  PAEQIEGCI
PKERS  235
```

ORDER INFORMATION

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.