



## Human Fibroblast Growth Factor 19 (FGF19). Recombinant

### Product Information

Code	00148-01-100
Name	Human FGF19, rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	Mature form
Protein ID	NP_005108
Gene ID	NM_005117
MW	22 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel PBS lyophilized
Formulation	Form without preservatives
Carry	free
Storage	-20 °C ~ -70° C
Reconstitution	500 µl
Application	ELISA

### Description

A DNA sequence encoding the mature form of human FGF19 (Pro24-Lys216) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

### Formulation

Lyophilized 100 µg of human FGF19 in 100 µl of PBS, Carry free.

### Reconstitution & Storage

Add 500 µl deionized water to the vial to prepare a working stock solution at 200 µg/mL. Allow to set at least 30 minutes at 4 °C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

**Sequence:** mature form of human FGF19 (Pro<sup>24</sup>-Lys<sup>216</sup>)

ORDER INFORMATION  
AVISCERA BIOSCIENCE INC.  
2348 Walsh Ave. Suite C  
Santa Clara, CA 95051  
USA  
Tel: (408) 982 0300  
Fax: (408) 982 0301  
Email:  
[Info@AvisceraBioscience.com](mailto:Info@AvisceraBioscience.com)  
[www.AvisceraBioscience.com](http://www.AvisceraBioscience.com)

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.