



# AVISCERA BIOSCIENCE

## Anti-Human Fetuin A IgG

### Product Information

Code	A00173-01-100
Name	Anti Human Fetuin A Antibody
Clone No.	N/A
Lot No.	
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	Human Fetuin A rec.
Ab Type	IgG
Purification	Protein A affinity
Formulation	lyophilized Form without preservatives
Carry	free
Storage	-20~ -70 ° C
Specificity	Human , Rat
Reconstitution	PBS, 100 µl
Application	IHC, ELISA

### PREPARATION

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant human Fetuin A. That IgG was purified by Protein A affinity

### FORMULATION

100 µg of Rabbit Anti Human Fetuin A purified IgG in 100 µl of PBS without preservatives was lyophilized.

### RECONSTITUTION

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for up six months without detectable loss of activity.

### STORAGE

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ~ -70 ° C for 12 months. **Avoid repeated freeze-thaw cycles.**

### SPECIFICITY

This antibody has been selected for its ability to recognize human Fetuin A in indirect ELISA and immunohistochemistry.

### APPLICATIONS

**Indirect ELISA** - This antibody can be used at 1:4000 (0.25 µg/ml) to detect recombinant Human Fetuin A on indirect ELISA.

**Immunohistochemistry**- That antibody can be used at 2-4µg/ml to detect Fetuin A in human hepatocellular carcinoma or rat liver tissue slides (ABC).

*Optimal dilutions should be determined by each laboratory for each application.*

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

### ORDER INFORMATION

AVISCERA BIOSCIENCE, INC  
2348 Walsh Ave., Suite C  
Santa Clara, CA 95051  
Tel: (408) 982 0300  
Info@Aviscerabioscience.com