



# AVISCERA BIOSCIENCE

## ANTI-BOVINE FETUIN A IGG

### Product Information

Code	A00173-07-100
Name	Bovine Fetuin A IgG
Clone No.	Polyclonal
Lot No.	
Size	100 µg
Species	Bovine
Host	Rabbit
Immunogen	Bovine Fetuin A rec.
Ab Type	IgG
Purification	Protein A
Formulation	lyophilized Form without preservatives free
Carry	
Storage	-20 ° C
Specificity	Bovine
Reconstitution	PBS, 100 µl
Application	IHC  ELISA

### Preparation

This antibody was produced from a rabbit immunized with purified bovine fetuin A. That was purified by Protein A affinity.

### Formulation

100 µg of rabbit anti-bovine Fetuin A IgG in PBS without preservatives was lyophilized.

### Reconstitution

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl.

Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted ( by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for up six months without detectable loss of activity.

### Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months.

**Avoid repeated freeze-thaw cycles.**

### Specificity

This antibody has been selected for its ability to recognize Bovine Fetuin A in direct ELISA.

### Applications

**Indirect ELISA** - This antibody can be used at 1:8000 (0.125µg/ml) with the appropriate secondary reagents to detect Bovine Fetuin A.

*Optimal dilutions should be determined by each laboratory for each application.*

AVISCERA BIOSCIENCE, INC  
2348 WALSH AVE., SUITE C  
SANTA CLARA, CA 95051  
USA

TEL: (408) 982 0300

Sales@Aviscerabioscience.com

www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

**AVISCERA BIOSCIENCE**