



AVISCERA BIOSCIENCE

Mouse Galectin-3 Recombinant

Product Information

Code	00199-03-100
Name	Galectin-3 (Mouse), rec.
Lot No.	
Size	100 µg
Species	Mouse
Sequence	Full length
Protein ID	NP_001139425
Gene ID	NM_001145953
MW	29 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS- PAGE gel PBS lyophilized
Formulation	Form without preservatives
Carry	free
Storage	-20 ° C ~ -70° C
Reconstituti on	500 µl
Application	ELISA

Description

A DNA sequence encoding the mouse galectin 3 (Met¹-Ile²⁶⁴) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 100 µg mouse Galectin-3 in 100 µl of PBS, Carry free.

Reconstitution & Storage

Add 500 µl PBS to the vial to prepare a working stock solution at 200 µg/mL . Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at - 20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: mouse Galectin-3 (Met1-Ile250)

AVISCERA BIOSCIENCE

2348 Walsh Ave., Suite C

Santa Clara, CA 95051

USA

Tel: (408) 982 0300

Fax: (408) 982 0301

Info@AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.