



# AVISCERA BIOSCIENCE

## Human Hypoxia and Akt induced Stem cell Factor (HASF) /C3ORF58 Recombinant

### Product Information

Code	00400-01-100
Name	Human HASF/C3ORF58 Rec
Lot No.	
Size	100 µg
Species	Human
Sequence	Leu36-Arg430
Protein ID	NP_775823.1
Gene ID	NM_173552.3
MW	48.9 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel
Formulation	PBS lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C ~ -70° C
Reconstitution	100 µl

### Description

A DNA sequence encoding the human HASF / C3ORF58 (Leu<sup>36</sup>-Arg<sup>430</sup>) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

### Formulation

Lyophilized 100 µg human HASF/C3ORF58 (36-430) in 100 µl of PBS, Carry free.

### Reconstitution & Storage

Add 100 µl PBS to the vial to prepare a working stock solution at 1000 µg/mL. Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months).

Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

**Sequence:** human HASF / C3ORF58 (Leu<sup>36</sup>-Arg<sup>430</sup>)

AVISCERA BIOSCIENCE

2348 Walsh Ave, Suite C

Santa Clara, CA 95051

USA

Tel: (408) 982 0300

Fax: (408) 982 0301

Info@AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.