



Product Information

Code	00539-01-100
Name	Human HMGB1 Rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	Pro9_Glu215
Protein ID	NP_002119
Gene ID	NM_002128
MW	27.9 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel PBS lyophilized form without preservatives
Formulation	
Carry	Free
Storage	-20 °C ~ -70° C
Reconstitution	100 µl
Application	Cell biology

Human High Mobility Group Protein B1 (HMGB1) Recombinant

Description

A DNA sequence encoding the mature form of human HMGB1 (Pro⁹-Glu²¹⁵) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 100 µg of human HMGB1 in 100 µl of PBS Carry free.

Reconstitution & Storage

Add 100 µl deionized water to the vial to prepare a working stock solution at 1000 µg/mL. Allow to set at least 30 minutes at 4 °C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: Human HMGB1 (Pro⁹-Glu²¹⁵)

9PR

```
GKMSSYAFFV QTCREEHKKK HPDASVNFSE FSKKCSERWK TMSAKEKGGF
EDMAKADKAR YEREMKTYIP PKGETKKKFK DPNAPKRPPS AFFLFCSEYR
PKIKGEHPGL SIGDVAKKLG EMWNTAADD KQPYEKKA AK LKEYEKDIA
AYRAKGKPPDA AKKGVVKA EK SKKKKEEEED EEEDEEEEEE EDEEDEDEEE
DDDDE 215
```