



AVISCERA BIOSCIENCE

Anti Irisin (Human, Mouse, Rat) Monoclonal Antibody

Product Information

Code	A00170-11-100
Name	Anti Irisin (H,M,R) MAB
Clone No.	
Lot No.	
Size	100 µg
Species	Human, Mouse, Rat
Host	Mouse
Immunogen	Irisin (Human, Mouse, Rat)
Ab Type	IgG
Purification	Protein G
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human,
Reconstitution	N/A
Application	IHC ELISA

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified recombinant Irisin (Human, Mouse, Rat). That antibody was purified by Protein G affinity.

Formulation

100 µg of Anti Irisin (Human, Mouse, Rat) Monoclonal Antibody in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (1000µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for six months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize Irisin (Human, Mouse, Rat) in direct ELISAs.

Applications

Indirect ELISA - This antibody can be used at 0.04 µg/ml to detect Irisin (Human, Mouse, Rat) on indirectly ELISA.

Immunohistochemistry- That Antibody can be used at 4-6 µg/ml with the appropriate secondary antibody to detect Irisin in human and rat skeletal muscle and heart tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE, INC
2348 Walsh Ave., Suite C
Santa Clara, CA 95051
USA

Tel: (408) 982 0300

Info@AvisceraBioscience.com

www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.