



AVISCERA BIOSCIENCE

Rabbit Anti-Human S100A6 IgG

Product Information

Code	A00316-01-100
Name	S100A6 (H) PAb,
Clone No.	Polyclonal
Lot No.	
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	S100A6 (H), rec.
Ab Type	Rabbit IgG
Purification	Protein A
Formulation	lyophilized Form without preservatives free
Carry	
Storage	-20 ° C
Specificity	H
Reconstitution	PBS, 100 µl
Application	IHC ELISA

Preparation

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant human S100A6 that was purified by Protein A.

Formulation

100 µg of rabbit anti-human S100A6 IgG in 100 µl PBS without preservatives was lyophilized.

Reconstitution

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize human S100A6 in direct ELISAs and western blots as well as immunohistochemistry.

Applications

Direct ELISA - This antibody can be used at 1:8000 (0.125 µg/ml) with the appropriate secondary reagents to detect human S100A6.

Immunohistochemistry-That Antibody can be used at 1: 500~ 1000 (1-2 µg/ml) with the appropriate secondary antibody to detect S100A6 in human lung cancer tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.