



AVISCERA BIOSCIENCE

Human Stathmin/Metablastin recombinant

Product Information

Code	00593-01-100
Name	Stathmin/Metabolastin (Human), rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	Mature form
Protein ID	NP_981946
Gene ID	NM_203401
MW	21.1 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel
Formulation	Tris buffer lyophilized Form without preservatives
Carry	free
Storage	-20 ° C ~ -70° C
Reconstitution	200 µl
Application	ELISA Cell biology

Description

A DNA sequence encoding the human Stathmin/Metabolastin (Ala2-Asp149) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 100 µg human Stathmin/Metabolastin in 100 µl of PBS pH7.4, Carry free.

Reconstitution & Storage

Add 200 µl deionized water to the vial to prepare a working stock solution at 500 µg/mL. Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months).

Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: human Stathmin/Metabolastin (Ala2-Asp149)

AVISCERA BIOSCIENCE, INC
2348 Walsh Ave., Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Info@AvisceraBioscience.com
WWW.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.