



Anti Human Vitamin D Binding Protein (VDBP) Monoclonal Antibody

Product Information

Code	A00678-01-100
Name	Anti Human VDBP Monoclonal Antibody
Clone No.	2B7B7B10
Lot No.	
Size	100 µg
Species	Human
Host	Mouse
Immunogen	Natural Human VDBP
Ab Type	IgG
Purification	Protein G
Formulation	Lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	100 µl
Application	ELISA , IHC

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified human plasma-derived VDBP. That antibody was purified by Protein G affinity.

Formulation

100 µg of Mouse anti human VDBP Monoclonal Antibody in 100 µl of PBS lyophilized form.

Reconstitution and Storage

Add 100 µl deionized water to the vial to prepare a antibody stocking solution (100µg/ml). Stores it at 4°C for a few days. For long term storage, the reconstituted antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for 12 months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize nature Human VDBP on indirect ELISAs.

Applications

Indirect ELISA - This antibody can be used at 1:8000 (0.125 µg/ml) to detect the natural human VDBP on indirectly ELISA.

Immunohistochemistry-That Antibody can be used at 1: 250-500 (2~4 µg /mL) with the appropriate secondary antibody to detect VDBP in the paraffin embedded human liver tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Info@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.