



# AVISCERA BIOSCIENCE

## Human Soluble VSTM2L (25-166) Recombinant

### Product Information

Code	00820-03-100
Name	Soluble VSTM2L (25-166) (Human) rec.
Lot No.	
Size	100 µg
Species	Human
Sequence	Mature form
Protein ID	Q96N03
Gene ID	128434
MW	19.83 KDa
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel PBS lyophilized
Formulation	Form without preservatives
Carry	free
Storage	-20 ° C ~ -70° C
Reconstitution	500 µl
Application	ELISA

### Description

A DNA sequence encoding the mature form of human Soluble VSTM2L (25-166) (Thr<sup>25</sup>-Gly<sup>166</sup>) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

### Formulation

Lyophilized 100 µg human Soluble VSTM2L (25-166) in 100 µl of PBS. Carry free.

### Reconstitution & Storage

Add 500 µl deionized water to the vial to prepare a working stock solution at 200 µg/mL. Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months). Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

**Sequence:** mature form of human Soluble VSTM2L (25-166)  
(Thr<sup>25</sup>-Gly<sup>166</sup>)

### AVISCERA BIOSCIENCE

2348 Walsh Ave., Suite C  
Santa Clara, CA 95051  
USA

Tel: (408) 982 0300

Email:

Info@Aviscerabioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.