



AVISCERA BIOSCIENCE

Human Asprosin / Fibrillin-1 (2732-2871) Recombinant

Product Information

Code	00229-01-50
Name	Human Asprosin Rec
Lot No.	
Size	50 µg
Species	Human
Sequence	Ser2732-2871His
Protein ID	NP_000129.3
Gene ID	NM_000138.4
MW	19.7 KD
Tag	His tag on N terminal
Source	E. Coli
Purity	>95% in SDS-PAGE gel
Formulation	PBS lyophilized Form without preservatives
Carry	Free
Storage	-20 ° C ~ -70° C
Reconstitution	100 µl of TBS, pH 8.0 , 9.0

Description

A DNA sequence encoding the human Asprosin /Fibrillin-1 (Ser²⁷³²-His²⁸⁷¹) with 6 His tag on the N-Terminus was expressed in *E. Coli*. This protein was purified by Ni-NTA column.

Formulation

Lyophilized 50 µg human Asprosin/ Fibrillin-1 (2732-2871) in 50 µl of PBS, Carry free.

Reconstitution & Storage

Add 100 µl TBS, pH 8.0 ~ 9.0 to the vial to prepare a working stock solution at 500 µg/mL. Allow to set at least 30 minutes at 4 ° C, mix well.

Store lyophilized protein at -20 °C or -70 °C. Lyophilized protein is stable for up to 6 months from date of receipt at -20 °C to -70 °C. Upon reconstitution, this protein can be stored at -20 °C for a few weeks or at -70 °C in a manual defrost freezer for long term storage (six months).

Aliquot reconstituted protein to avoid repeated freezing / thawing cycles.

Sequence: human Asprosin / Fibrillin-1 (Ser²⁷³²-His²⁸⁷¹)

AVISCERA BIOSCIENCE

2348 Walsh Ave, Suite C

Santa Clara, CA 95051

USA

Tel: (408) 982 0300

Fax: (408) 982 0301

Info@AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.