



AVISCERA BIOSCIENCE

Anti-Betatrophin/Lipasin/ANGPTL8 (Mouse) IgG Biotinylated

Product Information

Code	A00528-03-50B
Name	Betatrophin/Lipasin (Mouse) IgG
Clone No.	Polyclonal
Lot No.	
Size	50 µg
Species	Mouse
Host	Rabbit
Immunogen	Betatrophin (Mouse) rec.
Ab Type	IgG
Purification	Protein A
Formulation	lyophilized Form without preservatives free
Carry	
Storage	-20 ° C
Specificity	Mouse
Reconstitution	PBS, 50 µl
Application	IHC ELISA

Preparation

This antibody was produced from a rabbit immunized with purified recombinant Mouse Betatrophin/Lipasin/ANGPTL8. That was purified by Protein A affinity and conjugated with water soluble biotin.

Formulation

50 µg of rabbit anti-mouse Betatrophin/Lipasin IgG Biotinylated in 100 µl PBS contains 0.05% BSA without preservatives was lyophilized.

Reconstitution

Add 50 µl of PBS to the vial to prepare antibody stock solution.

Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C in a **manual defrost freezer** for up six months without detectable loss of activity.

Storage

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ° C for six months. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for its ability to recognize mouse betatrophin on indirect ELISAs and immunohistochemistry.

Applications

Indirect ELISA - This antibody can be used at 0.25 ~0.125 µg/ml with the appropriate secondary reagents to detect recombinant mouse betatrophin.

Immunohistochemistry-That Antibody can be used at 2-4 µg/ml with the appropriate secondary antibody to detect betatrophin/Lipasin in the paraffin embedded mouse liver and adipose tissues (ABC).

Optimal dilutions should be determined by each laboratory for each application.

AVISCERA BIOSCIENCE INC.
2348 Walsh Ave. Suite C
Santa Clara, CA 95051
USA
Tel: (408) 982 0300
Fax: (408) 982 0301
Email:
Sales@AvisceraBioscience.com
www.AvisceraBioscience.com

THIS PRODUCT IS FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.